N | Ryggstøtte, med OUTLAST®

1 | Levering

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GB | Shoulder support, with OUTLAST®

1 | Delivery

1 x Support belt for correcting posture

Read the instructions for use completely and carefully!

2 | Intended purpose

The support belt is designed to activate the abdominal and back muscles through external pressure (compression). At the same time, the shoulders are pulled back through the shoulder straps. The bandage straightens the upper body and relieves the lumbar spine. The stabilization and relief of the lumbar spine and the surrounding muscles relaxes and thus relieves pain. Supports the abdomen and back and promotes an upright posture

- Discomfort of the intervertebral discs
- Pain due to overloading
- Bad posture
- Lumbar vertebra syndrome
- Curvature of the spine

For private use

3 | Contraindications / Not suitable

In the case of the following medical conditions, it is only recommended to wear this article after consultation with your doctor:

- Skin diseases/injuries in the affected part of the body, especially in inflammatory symptoms such as swelling, redness and overheating
- Allergic reactions to the material
- Sensory and movement disorders, e.g. in the case of diabetes
- Open wounds or internal injuries
- For small children or persons who can not be noticeable when the product sits uncomfortably

4 | Notes on safe use of the product

Side effects affecting the whole organism are not yet known. Proper application is required. In case of improper use, product liability is excluded. The comfortable hook and loop fastener closes the straight holder by simply laying it on. There is no defect when small fibers are pulled from the surface material of the bandage when opening.

- If the bandage is too tight or uncomfortable, remove it immediately!
- Used by the times by private users.
- If your symptoms worsen or you notice unusual changes, consult your doctor immediately.

Do not use: while sleeping.

5 | Application

Fig. 1: Hold the support belt with the hook and loop closure (inside) in the right hand and guide it around the back. Pull the two straps sideways away from the body. Fig. 2: First fold in the left side and then the right one. The first closure allows a simplified system. Fig. 3: Then close the hook and loop fastener anywhere on the bandage. Lay the two shoulder straps over the shoulder, leading them below the axle and crossing them in the back area. Pull the straps forward on both sides and fix them at any point on the bandage. The stronger the pull, the stronger the support force.

6 | Cleaning

Cleaning is done by household laundry, according to care instructions. Close hook and loop fasteners during washing to avoid damaging other items of laundry.

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